



## Parent Instructions for Use and Care

Your child has received the InfantEar™ non-surgical ear molding system. Ear molding is not painful for your baby. In fact, they do not realize that the device is on.

During the 2 to 3 week course of treatment there are several things that can enhance your child's results:

Try to keep the beanie on as much as you can. It helps to protect the device from the child grabbing it or it potentially getting banged or knocked off.

The device should stay on for two-three weeks. You may notice the outer edges lifting slightly. That is ok, as long as main device remains intact. Your physician may provide you with an extra adhesive that you can apply to the edges if necessary.

Avoid getting the device wet - as this could loosen the adhesive and lead to the device coming off prematurely, or become reddened/irritated similar to a diaper rash.

Bathing - you can wash around the InfantEar device, but avoid any water dripping into the device.

Feeding – Breast or bottle feed normally, but be careful to avoid allowing formula, breastmilk or spit up to drip down into the device. If this happens, contact your physician to schedule an appointment to reapply the device as soon as possible.

Redness and/or Odor: You might notice some odor coming from the device. This is normal, as old skin cells are shed. However, if you see significant redness or seeping, please take a picture and contact your physician for a follow-up visit. This redness and/or odor is transient and will resolve soon after the device is removed.

For additional information on the InfantEar™ Non-surgical Ear Molding System please visit our website at [www.infantear.com](http://www.infantear.com)

---

TalexMedical, LLC 5 Great Valley Parkway Suite 328, Malvern PA 19335 1(888)EAR-2221

2-6-19 Rev. D